




ラルル

# Lulle

Heat & Cold Therapy Pack


Technology from **JAPAN**



<https://www.thaiolba.com/>



[info@thaiolba.com](mailto:info@thaiolba.com)



AD License No. XXXX/XXXX

BR-LUL2026EN-Rev01

# Sunflower Bean-Shaped Hot & Cool Pack



**Sherbet  
Gel**

Designed to fit  
comfortably on the  
chest and lymph nodes

Maintains optimal temperature

Fits snugly to the body

Gentle warmth to the body and behind the ears  
for relaxation and better sleep

Doubles as an eye mask



Insert the gel pack into the cotton cover



Notice the warning on the label and accompanying documents before use.



The cover: purple and gray

## How to prepare:

**Cold:** Place the gel pack in a freezer for at least 3 hours.

**Hot:** Heat the gel pack in a 500-600 watt microwave for 30 sec.

# Balloon Flower

## Intimate Hot & Cool Pack



Intimate  
Hot & Cool Pack







### How to prepare:

**Cold:** Place the gel pack in a freezer for at least 3 hours.

**Hot:** Heat the gel pack in a 500-600 watt microwave for 20 sec.

Insert the gel pack into the cotton cover



Gentle warmth  
to the uterine  
area

Relieve  
menstrual  
cramps

Reduce pain  
from postpartum  
vaginal wounds or  
hemorrhoids



The cover : purple and gray

Notice the warning on the label and accompanying documents before use.



# ラルル Lulle

Listing notification No. 69-2-3-2-0000628

## Main Components:

- **Film:** Nylon, polyethylene
- **Gel Contents:** Water, CMC, anti-freezing agent, control agent, preservative, colorant
- **Cover:** Cotton 100%



## Lulle Sunflower

Product Model	Weight (g)	Size (mm)	Colors
Lul-002	120	W90×H185	Purple Gray



## Lulle Ballon Flower

Product Model	Weight (g)	Size (mm)	Colors
Lul-003	90	W70×H180	Purple Gray

Importer: **THAI OLBA HEALTHCARE CO., Ltd.**  
131, Moo 9, Innovation Cluster 1, Room INC1-114,  
Khlong Nueng, Khlong Luang, Pathum Thani 12120 Thailand  
Tel (+66) 94-065-1409 , (+66) 94 059 9485  
<https://www.thaiolba.com/>

Manufacturer: **Mie Chemical Industry Co., Ltd.**  
255-1 Okuchi-chou Matsusaka-city Mie 515-0001 Japan  
TEL (0598) 51-2361 (代) FAX (0598) 51-1143  
<http://www.miekagaku.co.jp>

How to use



**THAI OLBA**